


# Wound Care Product Catalog

SOLUTIONS FOR  
HEALTH CARE  
PROVIDERS



# WOUND CARE

## INFECTION MANAGEMENT

Cutimed® Sorbact® – Hydrophobic Microbe Binding Antimicrobial Dressings.....	4
– Sorbact® Dressing Pad	
– Sorbact® Swab	
– Sorbact® Round Swab	
– Sorbact® Wound Contact Layer	
– Sorbact® Ribbon Gauze	
– Sorbact® Gel	
Cutimed® Sorbact® Hydroactive B .....	6
Cutimed® Siltec® Sorbact® .....	6
– Siltec® Sorbact® Sacral	
– Siltec® Sorbact® Oval	
Cutimed® Sorbion® Sorbact® .....	6

## EXUDATE MANAGEMENT

Cutimed® Sorbion® – Super-Absorbent Wound Dressing .....	7
– Sorbion® Sachet	
– Sorbion® Border	
– Sorbion® Border Sacrum	
– Sorbion® Sachet XL	
– Sorbion® Sachet S Drainage	
– Sorbion® Sachet Multi-Star	
– Sorbion® Sana Gentle	
– Sorbion® Sana Multi-Star	
Cutimed® Siltec® – Silicone Foam Dressing with Super-Absorbers .....	8
– Siltec® Silicone Foam Dressing	
– Siltec® Plus Foam Dressing	
– Siltec® L Foam Dressing	
– Siltec® B Foam Bordered Dressing	
– Siltec® B Oval	
– Siltec® Sacrum Foam Dressing	
– Siltec® 2D Heel Foam Dressing	
– Siltec® 3D Heel Foam Dressing	
Cutimed® Alginate – Calcium Alginate Dressing .....	9
Cutimed® Cavity – Absorbent Foam Dressing for Deep Wounds .....	9
Cutimed® Gel – Transparent Amorphous Hydrogel .....	10
Cutimed® Hydro L / Cutimed® Hydro B – Hydrocolloid Dressings.....	10
Cutimed® Hydrocontol – Hydropolymer Dressing .....	11

## COLLAGEN

Cutimed® Epiona® – Native Collagen Dressing .....	12
---	----

## POST-SURGICAL DRESSINGS

Leukomed® Sorbact® .....	13
Leukomed® Control – Transparent Post-Op Dressing with Absorbent Hydropolymer Gel Pad .....	13
Leukomed® T plus – Waterproof Adhesive Transparent Film Island Dressing .....	14
Leukomed® T – Waterproof Adhesive Transparent Film Dressing .....	14
Leukomed® – Non-Woven Adhesive Sterile Dressing with Absorbent Pad .....	15
Leukomed® IV – Waterproof Adhesive Transparent Sterile Film for IV.....	15
Leukomed® skin sensitive – For Wound Protection on Fragile or Sensitive Skin .....	16
Leukomed® T skin sensitive – Waterproof Adhesive Transparent Film Dressing .....	17
Leukomed® T plus skin sensitive – Waterproof Adhesive Transparent Film Island Dressing.....	17

## WIDE AREA FIXATION AND RETENTION

Hypafix® – Non-Woven Adhesive Fixation .....	18
Hypafix® skin sensitive – Soft and Secure Fixation for Fragile or Sensitive Skin .....	19
Hypafix® transparent – Waterproof Adhesive Transparent Film in Roll .....	19
Elastomull® – Non-Adhesive Gauze Fixation Bandage .....	20
Cover-Roll® Stretch – Stretch Adhesive Fixation Dressing .....	20

## FIRST AID DRESSINGS

Leukoplast® elastic – Stretchable Lightweight Fabric Dressing .....	21
---	----

## TUBULAR BANDAGES

Tensogrip® – Tubular Elastic Support Bandage.....	22
Tensoshape® – Elastic Tubular Bandage.....	22

## SKIN CARE

Cutimed® Protect Spray/Foam – Medical Skin Protection.....	23
Cutimed® ACUTE – Urea-Based Skin Care.....	23

## CONTACT LAYER

Cuticell® / Cuticell® Classic – Non-Adherent Dressing .....	24
Cuticell® Contact – Silicone Wound Contact Layer .....	24

## IRRIGATION PRODUCTS

SplashCap® – Wound Irrigation Shield .....	25
AbscessCap™ – Abscess & Deep Wound Irrigation Shield .....	25

NOTES .....	26
-------------	----

# CUTIMED® SORBACT®

## HYDROPHOBIC MICROBE BINDING ANTIMICROBIAL DRESSINGS



Other antimicrobial dressings on the market utilize alternative antimicrobials to kill bacteria or fungi, leaving behind bacterial cell debris and endotoxins. Cutimed® Sorbact® utilizes a physical mode of action versus a chemical mode of action. The hydrophobic interaction safely removes bacteria without cell death or negatively impacting the healing process.

### INDICATIONS




For all unclean, colonized and infected exudating wounds

- Chronic wounds, such as venous ulcers, diabetic ulcers and pressure ulcers
- Post-surgical wounds and dehisced wounds
- Traumatic wounds
- Wounds after excision of fistulas and abscesses

### FEATURES AND BENEFITS

- No healing delays due to bacterial death, cell debris or endotoxin release
- No known bacterial resistance has been described
- No cytotoxicity
- No known risk of allergies
- Sterile
- Not made with natural rubber latex

The dressing surface of Cutimed® Sorbact® is coated with a hydrophobic fatty acid derivate which attracts wound pathogens.. When the dressing is moist, wound bacteria and fungi are naturally drawn to the hydrophobic surface. They become irreversibly bound, unable to move back to the wound bed. The bound pathogens are removed with each dressing change, which helps clean contaminated, colonized and infected wounds by reducing the overall bacterial load.



Wound Depth	Wound Phase	Exudate
☐ Superficial	■ Necrotic	☒ Low
☑ Deep	■ Infected	☒ Low to moderate
☑ Superficial + deep	■ Sloughy	☒ Moderate to high
	■ Granulation	☒ High to Excessive
	■ Epithelialization	


### CUTIMED® SORBACT® DRESSING PAD



- Absorbent core
- For moderate to high exudate

Dimensions	Item #	UOM/Box
2.8 x 3.5 in. (7 x 8.9 cm)	7216111	5/box
2.8 x 3.5 in. (7 x 8.9 cm)	7216110	40/box
4 x 4 in. (10 x 10 cm)	7216213	5/box
4 x 4 in. (10 x 10 cm)	7216210	40/box
4 x 8 in. (10 x 20 cm)	7216307	20/box

### CUTIMED® SORBACT® WOUND CONTACT LAYER



- Impregnated acetate fabric
- Post-surgical, partial and full thickness wounds

Dimensions	Item #	UOM/Box
2 x 3 in. (5 x 7.6 cm)	7266200	10/box
4 x 4 in. (10 x 10 cm)	7266201	10/box
4 x 5 in. (10 x 12.7 cm)	7266202	10/box
4 x 8 in. (10 x 20 cm)	7266203	10/box
6 x 6 in. (15 x 15 cm)	7266204	10/box
8 x 8 in. (20 x 20 cm)	7266205	10/box

### CUTIMED® SORBACT® SWABS



- Flat, folded dressing
- For deep or superficial wounds

Dimensions	Item #	UOM/Box
1.6 x 2.4 in. (4 x 6.1 cm)	7216413	5/box
1.6 x 2.4 in. (4 x 6.1 cm)	7216412	40/box
2.8 x 3.5 in. (7 x 8.9 cm)	7216511	5/box
2.8 x 3.5 in. (7 x 8.9 cm)	7216510	40/box

### CUTIMED® SORBACT® RIBBON GAUZE



- Does not adhere to a moist wound bed
- For superficial or deep wounds / cavities

Dimensions	Item #	UOM/Box
0.8 x 19.7 in. (2 x 50 cm)	7216606	20/box
2 x 78.7 in. (5 x 200 cm)	7216705	10/box

### CUTIMED® SORBACT® ROUND SWABS



- Impregnated acetate round swabs
- For deep or superficial wounds

Description	Item #	UOM/Box
Round Swab	7216805	70 (14 x 5)

### CUTIMED® SORBACT® GEL DRESSING




- Used for slough or partly necrotic wounds
- Combines WCL with preservative-free hydrogel
- Helps promote autolytic debridement

Dimensions	Item #	UOM/Box
3 x 3 in. (7.6 x 7.6 cm)	7261112	10/box
3 x 6 in. (7.6 cm x 15 cm)	7261113	10/box



## CUTIMED® SORBACT® HYDROACTIVE B

- Hydrophobic microbe binding Sorbact® wound contact layer
- Hydropolymer gel matrix helps maintain a moist wound environment
- Minimizes the risk of maceration
- Ideal for absorption of thin exudate
- Fixation boarder keeps the dressing in place



Data of clinical observations (n = 100) indicated a noticeable change of phases of wound healing. Inflammatory phase is decreased and the proliferative phase is increased (granulation and epithelialization) after 10 days of Cutimed® Sorbact® Hydroactive treatment<sup>1</sup>.

<sup>1</sup> Hallern von, B., Absorbent DACC coated dressings, a new approach in treating wound infection: A multi-centre European study. Poster Presentation, EWMA 2012.



Dimensions	Item #	UOM/Box
2.8 x 3.3 in. (7 x 8.5 cm)	7993300	10/box
5.5 x 5.5 in. (14 x 14 cm)	7993302	10/box
7.5 x 7.5 in. (19 x 19 cm)	7993303	10/box
5.5 x 9.4 in. (14 x 24 cm)	7993304	10/box

## CUTIMED® SILTEC® SORBACT® MICROBE BINDING SILICONE FOAM DRESSING

- Hydrophobic microbe binding Sorbact® wound contact layer
- Super-absorbent hot-melt stripes ensure vertical absorption of exudate reducing the risk of maceration
- The moisture vapor transmission rate (MVTR) allows for dressing permeability while also maintaining a moisture balance



Dimensions	Item #	UOM/Box
3 x 3 in. (7.5 x 7.5 cm)	7325100	10/box
4 x 4 in. (10 x 10 cm)	7325116	10/box
5 x 5 in. (12.5 x 12.5 cm)	7325101	10/box
6 x 6 in. (15 x 15 cm)	7325102	10/box
7 x 7 in. (17.5 x 17.5 cm)	7325103	5/box
9 x 9 in. (22.5 x 22.5 cm)	7325104	5/box
Sacral 7 x 7 in. (17.5 x 17.5 cm)	7325105	5/box
Sacral 9 x 9 in. (22.5 x 22.5 cm)	7325106	5/box
Oval 3 x 4 in. (7 x 10 cm)	7325118	10/box

## CUTIMED® SORBION® SORBACT® HYDROPHOBIC MICROBE BINDING SUPER ABSORBENT DRESSING

- Binds and removes bacteria with a physical mode of action
- Absorbs large amount of exudate
- Retains exudate under compression, minimizing the risk of maceration
- Helps maintain a moist wound healing environment



Dimensions	Item #	UOM/Box
4 x 4 in. (10 x 10 cm)	7269808	10/box
4 x 8 in. (10 x 20 cm)	7269809	10/box
8 x 8 in. (20 x 20 cm)	7269810	10/box
8 x 12 in. (20 x 30 cm)	7269811	10/box

# CUTIMED® SORBION®

## SUPER-ABSORBENT WOUND DRESSINGS



Cutimed® Sorbion® dressings featuring **Hydration Response Technology™** offer a unique combination of mechanically modified cellulose fibers and gelling polymers, encased by a polypropylene outer layer. The interaction between the cellulose fibers and the gelling polymers helps regulate the wound environment for suitable wound bed preparation and tissue protection.



### SORBION® SACHET® S

Provides soft debridement while controlling excessive exudate

Dimensions	Item #	UOM/Box
3 x 3 in. (7.5 x 7.5 cm)	7323200	10/box
4 x 4 in. (10 x 10 cm)	7323206	10/box
4 x 8 in. (10 x 20 cm)	7323209	10/box
6 x 6 in. (15 x 15 cm)	7323212	10/box
8 x 8 in. (20 x 20 cm)	7323215	10/box
8 x 12 in. (20 x 30 cm)	7323218	10/box



### SORBION® BORDER

Combined silicone adhesive border for easier application and atraumatic dressing removal

Dimensions	Item #	UOM/Box
4 x 4 in. (10 x 10 cm)	7101000	10/box
6 x 6 in. (15 x 15 cm)	7101002	10/box
8 x 8 in. (20 x 20 cm)	7101003	10/box
4 x 9 in. (10 x 22.5 cm)	7101004	10/box
Sacral 7.6 x 8.6 in (19.5 x 22 cm)	7101008	10/box



### SORBION® SACHET® XL

Special size for larger wounds covering the entire lower leg

Dimensions	Item #	UOM/Box
10 x 18 in. (25 x 45 cm)	7324000	10/box

### INDICATIONS

- For moderate to high exuding wounds such as venous ulcers, diabetic ulcers or pressure ulcers

### FEATURES AND BENEFITS

- Superabsorber with high retention capacity
- Prevent surrounding skin from maceration
- Securely locks in exudate, even under compression
- Removes wound debris and locks in wound germs
- Maintains a moist wound environment by gel formation
- No binders or adhesives
- No risk of leakage, high tear strength
- Able to apply on both sides
- Odor reduction
- Sterile, atraumatic wound dressing
- Not made with natural rubber latex



### SORBION® SACHET® S DRAINAGE

Specially designed for catheters and drainage

Dimensions	Item #	UOM/Box
4 x 4 in. (10 x 10 cm)	7324100	10/box



### SORBION® SACHET® MULTI STAR

Fits perfectly to the anatomical shape of the body

Dimensions	Item #	UOM/Box
3 in. (8 cm)	7323701	10/box
5½ in. (12.5 cm)	7323703	10/box



### SORBION® SANA GENTLE

A 3-D structured contact layer made of polyethylene to protect granulation tissue

Dimensions	Item #	UOM/Box
3½ x 3½ in. (8.9 x 8.9 cm)	7323319	10/box
5 x 5 in. (12.5 x 12.5 cm)	7323322	10/box
9 x 5 in. (22.5 x 12.5 cm)	7323325	10/box
9 x 9 in. (22.5 x 22.5 cm)	7323314	10/box
13 x 9 in. (33 x 22.5 cm)	7323316	10/box



### SORBION® SANA MULTI STAR

Unique shape allows for flexible application of difficult-to-dress sites

Dimensions	Item #	UOM/Box
4½ in. (11 cm)	7323801	10/box
6½ in. (17 cm)	7323803	10/box

# CUTIMED® SILTEC®

## FOAM DRESSING WITH SILICONE WOUND CONTACT LAYER




Cutimed® Siltec® with its super -absorbent particles and silicone would contact layer, provide exudate management for moderate to high levels of exudate. It maintains a moist wound environment which supports moist wound healing. The silicone layer helps to minimize pain upon dressing removal.

### INDICATIONS

- For exuding wounds such as venous ulcers, diabetic ulcers or pressure ulcers

### FEATURES AND BENEFITS

- Reduces the risk of maceration
- Atraumatic dressing changes
- Visible saturation
- Strikethrough and leakage prevention
- Maintains a moist wound environment
- Sterile
- Not made with natural rubber latex
- For moderation to high levels of exudate



Product	Dimensions	Item #	UOM/Box
<b>Siltec Non-bordered</b> 	2 x 2 ½ in. (5 x 6 cm)	7328500	10/box
	4 x 4 in. (10 x 10 cm)	7328501	10/box
	4 x 8 in. (10 x 20 cm)	7328502	10/box
	6 x 6 in. (15 x 15 cm)	7328503	10/box
	8 x 8 in. (20 x 20 cm)	7328504	5/box
<b>Siltec PLUS Gentle tack</b> 	2 x 2 ½ in. (5 x 6 cm)	7328800	10/box
	4 x 4 in. (10 x 10 cm)	7328801	10/box
	4 x 8 in. (10 x 20 cm)	7328802	10/box
	6 x 6 in. (15 x 15 cm)	7328803	10/box
	8 x 8 in. (20 x 20 cm)	7328804	5/box
<b>Siltec L Thin, non-bordered</b> 	2 x 2 ½ in. (5 x 6 cm)	7328300	10/box
	4 x 4 in. (10 x 10 cm)	7328301	10/box
	6 x 6 in. (15 x 15 cm)	7328302	10/box
<b>Siltec B Adhesive Border</b> 	3 x 3 in. (7.5 x 7.5 cm)	7328400	10/box
	4 x 4 in. (10 x 10 cm)	7328405	10/box
	4 x 9 in. (10 x 22.5 cm)	7328406	10/box
	5 x 5 in. (12.5 x 12.5 cm)	7328401	10/box
	6 x 6 in. (15 x 15 cm)	7328402	10/box
	7 x 7 in. (17.5 x 17.5 cm)	7328403	5/box
	9 x 9 in. (22.5 x 22.5 cm)	7328404	5/box
<b>Siltec B Oval</b> 	Oval 3 x 4 in. (7.5 x 10 cm)	7328414	10/box
<b>Siltec Sacrum</b> 	7 x 7 in. (17.5 x 17.5 cm)	7328700	5/box
	9 x 9 in. (22.5 x 22.5 cm)	7328701	5/box
<b>Siltec Heel 2D</b> 	6 ½ x 9 ½ in. (16.5 x 24 cm)	7328805	5/box
<b>Siltec Heel 3D</b> 	6 ½ x 9 ½ in. (16.5 x 24 cm)	7328601	5/box



# CUTIMED® ALGINATE

## CALCIUM SODIUM ALGINATE DRESSING



Cutimed® Alginat is a sterile, non-medicated alginate dressing consisting of calcium and sodium alginate, with an integrity layer giving it high moisture stability for moderate to high exuding wounds.

### INDICATIONS

- For moderate to high exudate levels
- For venous ulcers, diabetic ulcers, pressure ulcers, post surgical wounds and donor sites.
- Suitable for use under compression
- For shallow and deep wounds

### FEATURES AND BENEFITS

- Highly absorbent and conformable
- Fast gelling and aids in autolytic debridement
- Maintains a moist wound environment
- Retains the exudate in a gel
- With sterile scissors it can be cut or it can be folded to fit a wound

Dimensions	Item #	UOM/Box
<b>Dressing</b>		
2 x 2" / 5 cm x 5 cm	7263408	10/box
4 x 4" / 10 cm x 10 cm	7263409	10/box
4 x 8" / 10 cm x 20 cm	7263410	10/box
<b>Rope</b>		
1 x 11.75" / 2.5 cm x 30 cm	7263411	5/box

# CUTIMED® CAVITY

## ABSORBENT FOAM DRESSING FOR DEEP WOUNDS



Cutimed® Cavity is a sterile wound dressing pad comprising a polyurethane hydrophilic open-cell foam.

### INDICATIONS

- For deep wounds such as venous ulcers, diabetic ulcers or pressure ulcers
- Wounds after excision of fistulas and abscesses
- Not indicated for dry wounds and cavity wounds with small openings

### FEATURES AND BENEFITS

- Non-adherent; allows for ease of dressing application and positioning in the wound and atraumatic dressing removal
- Easy to apply and remove
- Highly absorbent
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
2 x 2.75" / 5 cm x 7 cm	7262100	10/box
4 x 4" / 10 cm x 10 cm	7262101	10/box
6 x 6" / 15 cm x 15 cm	7262103	10/box

# CUTIMED® GEL

## TRANSPARENT AMORPHOUS HYDROGEL



**8 g**   **15 g**   **25 g**  
**0.3 oz**   **0.5 oz**   **0.8 oz**




Cutimed® Gel is a sterile clear, amorphous hydrogel with a pure formulation to reduce the risk of allergic reaction. Its superior moisture donation allows dry wounds to stimulate autolytic debridement by providing a moist wound environment which aids in the promotion of wound healing.

### INDICATIONS

- To stimulate and support autolytic wound debridement
- To rehydrate dry wounds and maintain a moist wound environment

### FEATURES AND BENEFITS

- No preservatives
- Icons on box to choose the right product
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
8 g / .3 oz. 	7261000	10/box
15 g / .5 oz. 	7261001	10/box
25 g / .8 oz. 	7261002	10/box

# CUTIMED® HYDRO L (LIGHT) / CUTIMED® HYDRO B

## HYDROCOLLOID DRESSINGS



Hydro B



Hydro L (Light)



Cutimed® Hydro L and Hydro B are sterile, self-adherent, hydrocolloid dressings. Cutimed® Hydro L is for superficial dry to light exuding wounds.

Hydro B has tapered edges, adhesive border, and a foam layer topped with a permeable film.

### INDICATIONS

- Hydro L (Light) – for dry to lightly exuding wounds
- Hydro B – for light to moderately exuding wounds

### FEATURES AND BENEFITS

- Facilitates and maintains a moist wound environment
- Forms gel in contact with exudate
- Vertical absorption of exudate
- Product composition allows removal of the gel with the dressing

Dimensions	Item #	UOM/Box
<b>Hydro L (Light)</b>		
2 x 10" / 5 cm x 25 cm	7263600	12/box
3 x 3" / 7.5 cm x 7.5 cm	7263601	10/box
4 x 4" / 10 cm x 10 cm	7263602	10/box
6 x 6" / 15 cm x 15 cm	7263603	5/box
<b>Hydro B</b>		
3 x 3" / 7.5 cm x 7.5 cm	7263500	5/box
4 x 4" / 10 cm x 10 cm	7263501	5/box
6 x 6" / 15 cm x 15 cm	7263502	5/box
Sacral - 6 x 7" / 15 cm x 17.5 cm	7263503	3/box

# CUTIMED<sup>®</sup> HYDROCONTROL

## HYDROPOLYMER DRESSING



Cutimed<sup>®</sup> HydroControl is an innovative combination of a hydropolymer gel matrix and a highly breathable polyurethane film. The interaction of these two components results in good fluid handling capacity and provides the ideal moist wound healing environment.

### INDICATIONS

- For use are on dry and low to moderately exuding wounds such as venous or arterial ulcers, diabetic foot ulcers or pressure ulcers

### FEATURES AND BENEFITS

- Promotes autolytic debridement and moist wound healing
- High wound stimulation potential
- For low to moderate exudate levels
- Optimum moisture balance

Dimensions	Item #	UOM/Box
1.75 x 1.75" / 4.5 cm x 4.5 cm	7323100	10/box
3 x 3" / 7.5 cm x 7.5 cm	7323101	10/box
4 x 4" / 10 cm x 10 cm	7323102	10/box
6 x 6" / 15 cm x 15 cm	7323103	10/box



# CUTIMED® EPIONA

FEATURING  
**3D MATRIX™ TECHNOLOGY**

**NATIVE COLLAGEN DRESSING -  
90% COLLAGEN AND 10% SODIUM ALGINATE**



Cutimed® Epiona is a native collagen, absorbent wound dressing designed to “jump start” the healing process in chronic and stalled wounds that have failed to meet their healing goals. A disrupted or damaged extracellular matrix (ECM) cannot support wound healing. The properties of a collagen wound dressing – approximating the skin’s native ECM structure – are acknowledged to be ideal in helping to activate healing.


- Better MMP binding by 50% over other collagen dressings
- Better protection of growth factors for an efficient wound-healing environment
- 7x more pronounced fibroblast cell growth supporting the formation of new granulation tissue



The open porous 3D matrix of Cutimed® Epiona acts as a scaffold to support the regeneration of human dermis that is absent from stalled, chronic wounds. A natural fenestrated structure of dermal extracellular matrix (ECM) collagen.<sup>1</sup>

- Reduces enzymatic degradation from inflammatory cells by binding MMPs
- Delivers structural support and proliferation of vital cells needed for wound repair<sup>2</sup>
- Creates an ECM-like scaffold encouraging cell proliferation and tissue growth

## Fibroblast Cell Growth



Fibroblast cell growth was studied in 14 days of laboratory testing on samples of Cutimed® Epiona, Promogran™ and Endoform™.

Cutimed® Epiona demonstrated more pronounced cell growth than the other collagens.



Cutimed® Epiona



Human Dermis

**Cutimed® Epiona’s native collagen provides a natural structure for enhanced cell growth.**

Dimensions	Item #	UOM/Box
2 x 2" / 5 cm x 5 cm	7322700	10/box
4 x 4" / 10 cm x 10 cm	7322701	10/box
8 x 8" / 20 cm x 20 cm	7322702	10/box

<sup>1</sup> Harding, K.; Kirsner, R.; Lee, D.; Mulder, G.; and Serena, T. International consensus: acellular matrices for the treatment of wounds. An expert working group review. Wounds International, 2010.

<sup>2</sup> Data on file.

# LEUKOMED® SORBACT®

**ANTIMICROBIAL WATERPROOF  
TRANSPARENT STERILE DRESSING  
WITH ABSORBENT PAD**



Leukomed Sorbact is a sterile, antimicrobial hydrophobic microbe binding post surgical island dressing. It utilizes a physical mode of action versus a chemical mode of action. It binds microbes to the sorbact layer and renders the bacteria inert

**INDICATIONS**

- Post-surgical incisions
- Post-operative dehisced wounds
- Lacerations, minor cuts, minor burns and abrasions

**FEATURES AND BENEFITS**

- Helps to minimize microbial colonization
- Effectively removes bacteria from the wound
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
2 x 3" / 5 cm x 7.2 cm	7619914	20/box
3 x 4" / 8 cm x 10 cm	7619915	20/box
3 x 6" / 8 cm x 15 cm	7619916	20/box
4 x 8" / 10 cm x 20 cm	7619917	20/box
4 x 10" / 10 cm x 25 cm	7619918	20/box
4 x 12" / 10 cm x 30 cm	7619919	20/box
4 x 14" / 10 cm x 35 cm	7619920	20/box

# LEUKOMED® CONTROL

**FULLY TRANSPARENT  
POST-OP DRESSING**



Leukomed® Control is a sterile transparent film dressing with an absorbent hydropolymer wound pad and a skin-friendly border. The transparency enables visual inspections of the wound without removing the dressing.

**INDICATIONS**

- Post-surgical incisions
- Post-operative dehisced wounds
- Acute wounds like minor cuts, minor burns and abrasions

**FEATURES AND BENEFITS**

- Creates optimal moisture level and adapts to skin's elasticity
- Innovative hydropolymer technology that manages exudate
- Fully transparent to help monitor for early detection of wound infection
- Convenient, all-in-one sterile dressing
- Good absorption capacity
- Skin friendly acrylic adhesive, secure yet easy to remove
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
2 x 2.75" / 5 cm x 7 cm	7323000	10/box
2.75 x 4" / 7 cm x 10 cm	7323001	10/box
3.125 x 6" / 8 cm x 15 cm	7323002	10/box
4 x 9.5" / 10 cm x 24 cm	7323003	5/box
4 x 13.75" / 10 cm x 35 cm	7323005	5/box



# LEUKOMED® T PLUS

## WATERPROOF ADHESIVE TRANSPARENT FILM ISLAND DRESSING



Leukomed® T plus is a sterile transparent film dressing with an absorbent wound pad for secreting wounds that is moisture permeable and waterproof.

### INDICATIONS

- Post-surgical and minimally invasive incisions
- For light to moderate wound drainage
- Lacerations, minor cuts, minor burns and abrasions

### FEATURES AND BENEFITS

- Easy to apply with the unique red stripe system
- Rounded edges prevent dressing from curling up
- Waterproof barrier to bacteria
- Permeable to water vapor and oxygen
- Convenient, ease of use, sterile dressing
- Sterile and not made with natural rubber latex

Dimensions	Item #	UOM/Box
2 x 3" / 5 cm x 7.2 cm	7238200	50/box
3 x 4" / 8 cm x 10 cm	7238201	50/box
3 x 6" / 8 cm x 15 cm	7238202	50/box
4 x 8" / 10 cm x 20 cm	7238214	50/box
4 x 10" / 10 cm x 25 cm	7238203	50/box
4 x 12" / 10 cm x 30 cm	7238204	50/box
4 x 14" / 10 cm x 35 cm	7238205	50/box

# LEUKOMED® T

## WATERPROOF ADHESIVE TRANSPARENT FILM DRESSING



Leukomed® T is a sterile transparent film dressing that is moisture permeable and waterproof. It is designed to cover wound sites, catheters or cannulas, while providing a barrier to external contaminants such as bacteria, blood and bodily fluids.

### INDICATIONS

- Superficial wounds
- Post-surgical and minimally invasive incisions
- Sterile fixation of wound dressings and catheters

### FEATURES AND BENEFITS

- Easy to apply with the unique red stripe system
- Rounded edges prevent dressing from curling up
- Waterproof and barrier to bacteria
- Permeable to water vapor and oxygen
- Completely transparent for easy inspection without removing dressing
- Sterile and single use dressing
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
2 x 3" / 5 cm x 7.2 cm	7238100	50/box
3 x 4" / 8 cm x 10 cm	7238101	50/box
4 x 4.75" / 10 cm x 12.5 cm	7238108	50/box
4.25 x 5.5" / 11 cm x 14 cm	7238102	50/box
6 x 8" / 15 cm x 20 cm	7238109	50/box
6 x 10" / 15 cm x 25 cm	7238110	50/box

# LEUKOMED®

## NON-WOVEN ADHESIVE STERILE DRESSING WITH ABSORBENT PAD



Leukomed® adhesive dressings are made from non-woven polyester fabric coated with a polyacrylate adhesive. This dressing is for short term use on the skin to cover minor and medium wounds after surgery and also minimally invasive surgeries. It can remain in place for one day to several days depending on the status of the wound and status of the skin.

### INDICATIONS

- For post-surgical and minimally invasive incisions
- For light to moderate drainage
- Lacerations, minor cuts, minor burns, and abrasions

### FEATURES AND BENEFITS

- Rounded edges prevent dressing from curling up
- Breathable
- For minor to moderate drainage
- Non-adherent wound pad doesn't stick to wound
- Conforms to the contours of the body
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
3 x 4" / 8 cm x 10 cm	7238001	50/box
3 x 6" / 8 cm x 15 cm	7238002	50/box
4 x 8" / 10 cm x 20 cm	7238003	50/box
4 x 10" / 10 cm x 25 cm	7238004	50/box
4 x 12" / 10 cm x 30 cm	7238005	50/box
4 x 14" / 10 cm x 35 cm	7238006	50/box

# LEUKOMED® IV

## WATERPROOF ADHESIVE TRANSPARENT STERILE FILM FOR I.V.



The Leukomed® IV transparent film dressing is for securing IV catheters and cannulas. This IV dressing allows visual IV site inspection and is bacteria proof, and is also moisture permeable.

### INDICATIONS

- Fixation of cannulae and peripheral venous catheters
- Puncture sites

### FEATURES AND BENEFITS

- Same advantages as Leukomed® T
- V-shaped design makes fixation easier (except child size)
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
1.75 x 1.75" / 4.5 cm x 4.5 cm	7239004	50/box
2 x 3" / 6 cm x 8 cm	7239000	50/box
2.75 x 3.5" / 7 cm x 9 cm	7239005	50/box
3.375 x 4.5" / 8.5 cm x 11.5 cm	7239003	50/box





# LEUKOMED® SKIN SENSITIVE

ABSORBENT DRESSING FOR  
WOUND PROTECTION ON  
FRAGILE OR SENSITIVE SKIN



Leukomed® skin sensitive absorbent dressings are made with a non-woven polyester fabric with a silicone adhesive boarder for ease of application. The skin sensitive silicone boarder allows for atraumatic dressing removal.

## INDICATIONS

- Post-surgical incisions
- Ideal for fragile and aging skin
- Acute wounds, minor cuts, minor burns, scrapes, and abrasions

## FEATURES AND BENEFITS

- Rounded edges prevent dressing from curling up
- Skin friendly with atraumatic removal
- Provides protection for closed surgical wounds that includes minimally invasive surgical wounds
- Absorbs blood and wound fluid
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
2 x 3" / 5 cm x 7.2 cm	7123904	5/box
3 x 4" / 8 cm x 10 cm	7123905	5/box
3 x 6" / 8 cm x 15 cm	7123906	5/box
4 x 10" / 10 cm x 25 cm	7123907	5/box

# LEUKOMED® T SKIN SENSITIVE

## WATERPROOF ADHESIVE TRANSPARENT FILM DRESSING



Leukomed® T skin sensitive is a transparent film dressing with a silicone adhesive boarder for secure wound protection and fixation on fragile skin. The skin sensitive silicone boarder allows for atraumatic dressing removal.

### INDICATIONS

- Post-surgical incisions
- Ideal for fragile and aging skin
- Acute wounds, minor cuts, minor burns, scrapes, and abrasions

### FEATURES AND BENEFITS

- Bacteria barrier and waterproof protection of post-operative and superficial traumatic wounds
- Sterile fixation of primary dressings and compresses
- Additional fixation of devices
- Protection of non-secreting wounds from friction and shear forces
- Sterile

Dimensions	Item #	UOM/Box
2 x 3" / 5 cm x 7.2 cm	7646400	5/box
3 x 4" / 8 cm x 10 cm	7646401	5/box
3 x 6" / 8 cm x 15 cm	7646402	5/box

# LEUKOMED® T PLUS SKIN SENSITIVE

## WATERPROOF ADHESIVE TRANSPARENT FILM ISLAND DRESSING



Leukomed® T plus skin sensitive transparent film island dressings with a silicone adhesive boarder for secure wound protection on fragile skin. The skin sensitive silicone border allows for atraumatic dressing removal.

### INDICATIONS

- Post-surgical incisions
- Acute wounds, minor cuts, minor burns, scrapes, and abrasions

### FEATURES AND BENEFITS

- Bacteria barrier and waterproof protection of surgical incisions
- Covers and protects minor cuts, minor lacerations and minor burns
- Absorbs blood and wound secretions
- Sterile

Dimensions	Item #	UOM/Box
2 x 3" / 5 cm x 7.2 cm	7617800	5/box
3 x 4" / 8 cm x 10 cm	7617801	5/box
3 x 6" / 8 cm x 15 cm	7617802	5/box
4 x 10" / 10 cm x 25 cm	7617803	5/box





# HYPAFIX®

## NON-WOVEN WIDE AREA FIXATION



Hypafix® is indicated for wide-area fixation and widely used in place of surgical tapes to secure and retain primary dressings on surgical or wound sites.

### INDICATIONS

- Fixation of primary dressing on small to large surgical or wound sites
- Ideal use for areas of body requiring flexibility such as elbows and knees
- Ideal base coverage for areas applying the patello-femoral taping technique

### FEATURES AND BENEFITS

- Soft, stretchable non-woven polyester conforms well to the body
- Secure for up to 7 days
- Breathable
- Skin friendly adhesive
- Easily cut to eliminate waste
- Not made with natural rubber latex

Dimensions	Item #	Local Code	UOM/Box
2" x 11 yd / 5 cm x 10 m	7144301	4209	24/box
4" x 11 yd / 10 cm x 10 m	7144302	4210	24/box
6" x 11 yd / 15 cm x 10 m	7144303	4211	24/box
2" x 2 yd / 5 cm x 1.8 m	7193921	4215	36/box
4" x 2 yd / 10 cm x 1.8 m	7193919	4216	12/box
6" x 2 yd / 15 cm x 1.8 m	7193920	4217	12/box



# HYPAFIX® SKIN SENSITIVE

**SOFT & SECURE FIXATION FOR  
FRAGILE OR SENSITIVE SKIN**



Hypafix® skin sensitive is designed especially for use on patients with fragile or compromised skin. The soft silicone adhesive helps to protect against adhesive-related injuries.

### INDICATIONS

- Ideal for aging, fragile skin
- For highly contoured parts of the body
- Fixation of dressings, catheters, cannulae, patches and ointment

### FEATURES AND BENEFITS

- Soft, silicone, skin-friendly adhesive
- Atraumatic removal
- Easy repositioning
- Not made with natural rubber latex
- Convenient roll format

Dimensions	Item #	UOM/Box
2" x 5.5 yd / 5 cm x 5 m	7996601	1/box
4" x 5.5 yd / 10 cm x 5 m	7996602	1/box
4" x 2.5 yd / 10 cm x 2 m	7996604	1/box

# HYPAFIX® TRANSPARENT

**WATERPROOF ADHESIVE  
TRANSPARENT FILM IN ROLL**



The innovative waterproof film of Hypafix® transparent is very thin and stretchable while maintaining superior strength and durability. It also feels soft and adapts well to the body contours like a second skin.

### INDICATIONS

- Skin protection against friction and shearing
- Fixation of dressings, catheters, cannulae, patches and ointment

### FEATURES AND BENEFITS

- Easy to apply with the unique red stripe system
- Waterproof
- Permeable to water vapour and oxygen
- Completely transparent for easy inspection without removing dressing
- Convenient roll format
- Not made with natural rubber latex

Dimensions	Item #	UOM/Box
2" x 11 yd / 5 cm x 10 m	7238700	12/box
4" x 11 yd / 10 cm x 10 m	7237801	12/box
6" x 11 yd / 15 cm x 10 m	7237802	12/box
4" x 2.1 yd / 10 cm x 2 m	7237803	12/box

# ELASTOMULL®

## ELASTIC GAUZE BANDAGE



Elastomull® offers a precision fit, even in areas where there is constant movement such as the elbow or knee.

### INDICATIONS

- Fixation of dressings

### FEATURES AND BENEFITS

- High density fabric for a soft, rich texture
- Performance and stretchability to help control bleeding and edema
- Precision fit in areas with constant movement
- Helps ensure no slippage once bandage is in place
- Not made with natural rubber latex

# COVER-ROLL®

## STRETCH

## CONFORMING ADHESIVE FIXATION DRESSING



Cover-Roll® Stretch retains the primary dressing on small to large wound sites and is ideal to use on areas of body requiring flexibility.

### INDICATIONS

- Fixation of primary dressing on small to large surgical or wound sites
- Ideal use for areas of body requiring flexibility such as elbows and knees
- Ideal base coverage for areas applying the patello-femoral taping technique

### FEATURES AND BENEFITS

- Soft, stretchable non-woven polyester material adapts well to body contours
- Provides light compression without causing constriction or congestion
- The polyacrylate adhesive holds dressing safely and reliably in place
- Not made with natural rubber latex

Dimensions	Item #	Local Code	UOM/Box
<b>Non-sterile</b>			
1" x 4.1 yd	7780200	02088000	24/box
2" x 4.1 yd	7780300	02089000	12/box
3" x 4.1 yd	7780400	02101000	12/box
4" x 4.1 yd	7780500	02102000	12/box
<b>Sterile</b>			
1" x 4.1 yd	7779700	02075001	12/box
2" x 4.1 yd	7779800	02076001	12/box
3" x 4.1 yd	7779900	02070001	12/box
4" x 4.1 yd	7780000	02071001	12/box

Dimensions	Item #	Local Code	UOM/Box
2" x 2 yd / 5 cm x 1.8 m	4554700	45547	1/box
4" x 2 yd / 10 cm x 1.8 m	4554800	45548	1/box
6" x 2 yd / 15 cm x 1.8 m	4554900	45549	1/box
8" x 2 yd / 20 cm x 1.8 m	4555000	45550	1/box
12" x 2 yd / 30 cm x 1.8 m	4555100	45551	1/box
2" x 10 yd / 5 cm x 9 m	4555200	45552	1/box
4" x 10 yd / 10 cm x 9 m	4555300	45553	1/box
6" x 10 yd / 15 cm x 9 m	4555400	45554	1/box
8" x 10 yd / 20 cm x 9 m	4555500	45555	1/box
12" x 10 yd / 30 cm x 9 m	4555600	45556	1/box

# LEUKOPLAST® ELASTIC

## STRETCHABLE LIGHTWEIGHT FABRIC DRESSING



Leukoplast® elastic fabric adhesive bandages are made to absorb liquid fast. Featuring an extra-large wound pad capable of absorbing 10 times its own weight in water, the 360° adhesive surrounding the island pad helps seal the afflicted area.

### INDICATIONS

- Minor acute wounds like cuts, lacerations, abrasions

### FEATURES AND BENEFITS

- Available in 13 different shapes and sizes
- Perfect for use on knees, elbows, chins, noses and knuckles
- Elastic fabric moves with the curves and contours of body
- Sterile
- Not made with natural rubber latex

Dimensions		Item #	UOM/Box
¾ x 3" strip		7645617	100/box
1 x 3" strip		7645606	100/box
1 x 3" strip		7645604	300/box
¾" round		7645616	100/box
1¼" round		7645613	100/box
1½ x 3" patch		7645612	100/box
2 x 3" patch		7645610	50/box
4 x 2¾" patch		7645605	50/box
1½ x 3" knuckle		7645607	100/box
2½ x 1½" toe shield		7645618	100/box
Sm finger tip		7645608	100/box
Lg finger tip		7645609	50/box
3 x 3" 4-wing		7645611	50/box

**FORMERLY KNOWN AS Coverlet®**



Description	Item #	UOM/Box
Cabinet (empty)	7253702	1 unit
Cabinet (filled)	7253703	1 box per drawer
Eye Occluder - Regular	7646600	30/box
Eye Occluder - Junior	7646700	30/box

# TENSOGRIP®\*

## TUBULAR ELASTIC SUPPORT BANDAGE



Tensogrip® is an elasticized, tubular support bandage, ideal for the effective management of strains, sprains and soft tissue injuries. It provides pressure and support over varying body contours.

### INDICATIONS

- Light support and compression
- Retention of wound dressings

### FEATURES AND BENEFITS

- Provides light support, light compression and dressing retention
- Horizontal stretch
- Easily cut to size
- Air-permeable

Dimensions	Item #	Local Code	UOM/Box
1.75 in x 11 yd	7151300	7579	1/box
2.5 in x 11 yd	7151400	7580	1/box
2.75 in x 11 yd	7151900	7581	1/box
3 in x 11 yd	7151500	7582	1/box
3.5 in x 11 yd	7152000	7583	1/box
4 in x 11 yd	7151600	7584	1/box
4.5 in x 11 yd	7152100	7585	1/box
7 in x 11 yd	7151700	7586	1/box
8.5 in x 11 yd	7152300	7587	1/box
2.5 in x 11 yd Full Leg	7151401	7580FL	1/box
2.75 in x 11 yd Full Leg	7151901	7581FL	1/box
3 in x 11 yd Full Leg	7151501	7582FL	1/box
3.5 in x 11 yd Full Leg	7152001	7583FL	1/box
4 in x 11 yd Full Leg	7151601	7584FL	1/box
4.5 x 11 yd Full Leg	7152101	7585FL	1/box

# TENSOSHAPÉ®\*

## TUBULAR ELASTIC BANDAGE



Tensoshape® is an elastic tubular bandage suitable for light support. Its anatomical design provides graduated compression, which promotes improved blood flow and lymphatic return. Both full leg and below-the-knee bandages are available in a variety of sizes.

### INDICATIONS

- Light support and compression
- Treatment of minor soft tissue injuries

### FEATURES AND BENEFITS

- Easy to apply
- Anatomical design
- Applies graduated compression

Size	Calf Size	Item #	Local Code	UOM/Box
<b>Below the Knee</b>				
Small (B/C)	12.5 x 14"	7208502	7589	10/case
Medium (C/D)	13.7 x 15.2"	7208501	7590	10/case
Large (D/E)	14.8 x 16.4"	7208500	7591	10/case
<b>Full Leg</b>				
Small (B/C)	12.5 x 14"	7208506	7592	10/case
Medium (C/D)	13.7 x 15.2"	7208505	7593	10/case
Large (D/E)	14.8 x 16.4"	7208504	7594	10/case
X-Large (D/E)	16 x 17.6"	7208503	7595	10/case

\* Caution: These products contains natural rubber latex which may cause allergic reactions.

# CUTIMED® PROTECT

## MEDICAL SKIN PROTECTION



Cutimed® PROTECT provides a long lasting protective barrier film that shields at risk skin from irritants, such as exudate or incontinence.

### FEATURES AND BENEFITS

- Protects wound margins and peri-wound skin
- Preserves and maintains intact skin
- Prevents maceration and peri-wound skin breakdown
- Promotes dressing adhesion and pain free dressing removal
- Available in foam applicator and spray

Dimensions	Item #	Local Code	UOM/Box
<b>Spray</b>			
28 ml / .95 oz.	7265300	7265300	1/box
28 ml / .95 oz.	7265301	7265301	12/box
<b>Foam Applicator</b>			
1 ml / .03 oz.	7265400	7265400	5/box
1 ml / .03 oz.	7265401	7265401	25/box
3 ml / 1.1 oz.	7265402	7265402	5/box
3 ml / 1.1 oz.	7265403	7265403	25/box

# CUTIMED® ACUTE

## UREA-BASED SKIN CARE



Cutimed® ACUTE is a cream mousse which brings valuable relief to dry-skin sufferers and simple, practical assistance to their healthcare professionals. This high quality range of urea-based products offers medical skin care and protection for all types of dry skin, particularly stressed or elderly skin.

### INDICATIONS

- Protection for stressed or elderly skin

### FEATURES AND BENEFITS

- Absorbs rapidly, penetrates deeply
- Rebalances skin moisture levels
- Re-establishes epidermal barrier function
- Foam mousse is gentle to apply
- Leaves skin smooth and supple
- Does not leave greasy residue
- Matches urea concentration to patient's needs

	Dimensions	Item #	Local Code	UOM/Box
5%	125 ml / 4.23 oz.	7264122	7264122	1/box
10%	125 ml / 4.23 oz.	7264123	7264123	1/box





# CUTICELL® CONTACT

## SILICONE WOUND CONTACT LAYER



Cuticell® Contact is a silicone wound contact layer that has a perforated transparent structure and is highly conformable. Its optimally perforated structure ensures even highly viscous exudate can easily pass through.

### INDICATIONS

For superficial exuding wounds such as

- First and second-degree burns
- Lacerations and abrasions
- Radiation injuries
- Donor and recipient skin graft sites

### FEATURES AND BENEFITS

- Low-adherent and soothing to minimize pain at dressing change
- Supports a moist environment to help healing
- Allows exudate to drain freely into an absorbent secondary dressing
- High transparency for easy wound inspection
- Reduced risk of maceration
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	Local Code	UOM/Box
2 x 3" / 5 cm x 7.5 cm	7268000	7268000	5/box
3 x 4" / 7.5 cm x 10 cm	7268001	7268001	5/box
4 x 7" / 10 cm x 18 cm	7268002	7268002	5/box
6 x 10" / 15 cm x 25 cm	7268003	7268003	5/box

# CUTICELL®

## NON-ADHERENT DRESSING IMPREGNATED WITH NON-MEDICATED OINTMENT



Cuticell® is a low-adherent primary wound contact layer made of smooth small meshed acetate fabric impregnated with non-medicated ointment basis. The result: a dressing that ensures the wound is not disturbed and helps maintain a moist environment that supports wound healing.

### INDICATIONS

For superficial exuding wounds such as

- First and second-degree burns
- Lacerations and abrasions
- Radiation injuries
- Donor and recipient skin graft sites

### FEATURES AND BENEFITS

- Low-adherent and soothing to minimize pain at dressing change
- Supports a moist environment to help healing
- Does not disturb wound healing
- Does not charge the wound with excess ointment
- Simple handling without release paper
- Sterile
- Not made with natural rubber latex

Dimensions	Item #	Local Code	UOM/Box
3 x 3" / 7.5 cm x 7.5 cm	7253900	7253900	10/box
3 x 3" / 7.5 cm x 7.5 cm	7253901	7253901	50/box
3 x 8" / 7.5 cm x 20 cm	7253902	7253902	10/box
3 x 8" / 7.5 cm x 20 cm	7253903	7253903	50/box

# SPLASHCAP®

## WOUND IRRIGATION SHIELD



SplashCap® is a twist-on wound irrigation shield that allows a clinician to deliver hundreds of milliliters of irrigation fluid directly to a wound in seconds.


### FEATURES AND BENEFITS

- Faster wound irrigation improved productivity and faster turnaround times
- Built-in nozzle avoids the risk of needle sticks from poking holes in bottles
- Built-in shield reduces splattering when using bottle caps without shields
- Proven pressure (PSI): 4-15 PSI
- Universal fitting - fits all major standard saline bottles
- High volume - delivers over 250 ml of irrigation fluid in seconds
- Kid friendly - no scary syringes

Item #	UOM
7446800	25 box/case
7446801	100 box/case

# ABSCCESSCAP™

## ABCESS & DEEP WOUND IRRIGATION SHIELD



AbscessCap™ is a twist on abscess irrigation shield that allows a clinician to easily irrigate wounds below the surface of the skin.

### FEATURES AND BENEFITS

- Protruding nozzle - allows penetration below the surface
  - one-handed probing
  - break up loculations
  - irrigate tunneling wounds
- Fits any standard saline bottle - delivers high volumes in a shorter amount of time
- Fits any standard syringe - versatility between preferred irrigation methods
- Built in splash shield permits safe visualization of wound target while reducing splatter
- Large nozzle opening creates a low pressure stream of fluid to safely flush the wound of pus and debris

Item #	UOM
7446806	50 box/case

# NOTES



AbscessCap™	25
Cover-Roll® Stretch	20
Cuticell® / Cuticell® Classic	24
Cuticell® Contact	24
Cutimed® Acute	23
Cutimed® Calcium Alginate	9
Cutimed® Cavity	9
Cutimed® Epiona	12
Cutimed® Gel	10
Cutimed® Hydro L / Hydro B	10
Cutimed® Hydrocontrol	11
Cutimed® PROTECT	23
Cutimed® Siltec®	8
Cutimed® Siltec® Sorbact®	6
Cutimed® Sorbact®	4-5
Cutimed® Sorbact® Hydroactive B	6
Cutimed® Sorbion®	7
Cutimed® Sorbion® Sorbact®	7
Elastomull®	20
Hypafix®	18
Hypafix® skin sensitive	19
Hypafix® transparent	19
Leukoplast® elastic	21
Leukomed®	15
Leukomed® Control	13
Leukomed® IV	15
Leukomed® skin sensitive	16
Leukomed® Sorbact®	13
Leukomed® T	14
Leukomed® T plus	14
Leukomed® T plus skin sensitive	17
Leukomed® T skin sensitive	17
SplashCap®	25
Tensogrip®	22
Tensoshape®	22