

CUTIMED® OFF-LOADER SELECT

Total Contact Cast System

THE
**GOLD
STANDARD**
IN DIABETIC FOOT
ULCER HEALING

**TAKE THE
LOAD OFF.
TURN THE
HEALING ON.**

- TRUE total contact cast
- Increased patient compliance
- Locks in toe and heel
- Customizable for most legs
- Cost-effective

Cutimed®

CUTIMED® OFF-LOADER SELECT

Evidence-based treatment guidelines agree the Total Contact Cast (TCC) System is the Gold Standard for off-loading, which is important for healing of diabetic foot ulcers¹. Essity, a leading global hygiene and health company, brings you the Cutimed® Off-Loader Select Total Contact Casting System in an easy to use kit.

¹Consensus recommendations, Snyder, Kirsner, Laverly, OWM 2010

Guidelines for the treatment of diabetic ulcers, Steed et al. Wound Rep Reg (2006) Wound Healing Society




Evidence Based Results

Total contact casting provides higher percentage of ulcers healed vs other methods of offloading.²

- 90% TCC
- 65% RCW (Removable Cast Walker)
- 58% Half Shoe

²Armstrong Diabetes Care 2001; 24:1019-22. [RCT]

Off-Loading - Comparing Healing Rates Against Methods



²Armstrong Diabetes Care 2001; 24:1019-22. [RCT]

Total Contact Cast System

Benefits of Using Cutimed® Off-Loader Select Total Contact Cast System

- True total contact
- Customizable and adaptable to structural foot deformities with no weight limit
- Locks the heel and toes in place to significantly reduce shear, stress and friction
- Works to significantly reduce repetitive injury¹
- Comes in complete, easy to fit kits that take less than 10 minutes to apply
- Helps achieve better patient compliance¹



¹Total contact cast wall load in patients with a plantar forefoot ulcer and diabetes, Begg, et al. Journal of Foot and Ankle Research (2016) 9:2 DOI 10.1186/s13047-015-0119-0

Comparison of Cutimed® Off-Loader Select to Other Therapies

	% Healed	Healing Time
Total Contact Cast ²	80-90%	42 days
Removable Boot ²	52-74%	38-50 days
Apligraf ^{®3}	56%	65 days
Dermagraft ^{®4}	50%	72 days
Regranex ⁵	50%	140 days

²Bus, S.A., Valkm G.D. The effectiveness of footwear and offloading interventions to prevent and heal foot ulcers and reduce plantar pressure in diabetes: a systematic review, Diabetes Metab Res Rev 2008; 24(Suppl 1): S162-S180.

³Zaykyanov, Larissa; Kirsner, Robert S., A review of a bi-layered living cell treatment (Apligraf[®]) in the treatment of venous leg ulcers and diabetic foot ulcers, Clinical Interventions in Aging 2007;2(1) 93-98.

⁴Genzkow G.D., Iwasaki S.D., Hershon K.S., et al. Use of Dermagraft, a cultured human dermis, to treat diabetic foot ulcers. Diabetes Care, 1996 Apr;19:350-4

⁵Fang, Robert C., Galiano, Robert D., A review of becaplermin gel in the treatment of diabetic neuropathic foot ulcers.

Why a kit?

- Time savings
- Standardize treatment
- Convenience
- May help reduce inventory requirements



KIT COMPONENTS



Cutimed® Cavity is a sterile dressing, made from hydrophilic absorbent polyurethane foam. It is soft, highly absorbent and conformable creating the ideal environment for moist wound healing.



Delta Terry-Net™ C cotton stockinette has excellent absorption and conformability. Reduces the risk of bunching in the dorsum of the foot.



Delta Terry-Net™ adhesive felt padding for added protection of the bony prominences of the tibia and malleoli.



Specialist® 100 cotton cast padding is soft and made of 100% cotton. It has excellent absorption properties and low loft needed for a total contact cast.



Delta-Lite® Conformable is a patented fiberglass cast tape that provides excellent conformability, greatly reducing the need to tuck and fold during application.



Perforated open cell foam to protect the toes and allow more breathability than nonperforated foams.



Cotton gauze used between the toes to reduce the risk of skin maceration.

Cutimed® Off-Loader Select - Synthetic

Ref.-No.	Size	Qty
7800901	Cutimed® Cavity	1 each
	2" x 2" Cotton Gauze	4 each
	3" Delta-Lite® Conformable	3 rolls
	4" Delta-Lite® Conformable	2 rolls
	3" Delta Terry-Net™ Stockinette	1 each
	4" Specialist® 100 Cotton Cast Padding	2 each
	Delta Terry-Net™ Adhesive Felt	1 each
	Open Cell Foam	1 each

Canvas Rocker Bottom Cast Shoe

Ref.-No.	Description	Qty
43101408	Canvas Rocker Bottom Cast Shoe - Small	1 each
43101505	Canvas Rocker Bottom Cast Shoe - Medium	1 each
43101602	Canvas Rocker Bottom Cast Shoe - Large	1 each
43101709	Canvas Rocker Bottom Cast Shoe - X-Large	1 each

Cutimed® Off-Loader Select - Plaster & Synthetic Combination

Ref.-No.	Size	Qty
7800900	Cutimed® Cavity	1 each
	2" x 2" Cotton Gauze	4 each
	4" Orthoflex® Elastic Plaster*	3 rolls
	4" Extra Fast Gypsona® S	2 rolls
	3" Delta-Lite® Conformable	1 each
	4" Delta-Lite® Conformable	2 each
	3" Delta Terry-Net™ Stockinette	1 each
	4" Specialist® Cotton Cast Padding	1 each
	Delta Terry-Net™ Adhesive Felt	1 each
	Open Cell Foam	1 each

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

*Caution: This product contains natural rubber latex which may cause allergic reactions.

Procedure Coding

CPT/APC Code Descriptor
29445/5102 Application of the ridge
Total Contact Leg Cast

Facility Reimbursement

CPT 29445 Links to APC Code 5102

Physician Reimbursement

CPT 29445

REIMBURSEMENT DISCLAIMER: The Centers for Medicare & Medicaid Services (CMS) uses Palmetto GBA as the Pricing, Data Analysis, and Coding (PDAC) Contractor. Their website, located at <https://www.dmeopdac.com>, maintains a searchable database for HCPCS codes. Essity, however, makes no representations as to the accuracy of the information contained within the PDAC database, nor any representations as to whether its products are reimbursable under any government sponsored healthcare program and/or private-insurance program. It is the provider's sole responsibility, in consultation with the insurer, if necessary, to determine and submit appropriate codes, charges, and modifiers for services rendered.



essity

Cutimed®,
an Essity brand



/BSNWoundCareUSA



@BSNWoundCareUSA



bsnmedical.us

BSN Medical Inc., an Essity company
5825 Carnegie Blvd., Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325
To order toll-free: Tel. (+1) 800 537 1063

60734 R6 ©2019 BSN Medical Inc. REV 08/19