



medi

## Compression Solutions for Venous Leg Ulcers

Treating the underlying condition of CVI

medi. I feel better.

# circaid® juxtalite®

The clinically proven  
compression wrap<sup>1</sup>




Easy to use


Measurable  
Compression

Patient  
Independence

## NEW EZ-on Panel\*



Can be used day  
and night for  
maximum therapy



Measurable  
compression for  
proper dosage




Instantly adjustable  
for enhanced fit with  
juxtaposed bands



Clima Fresh  
prevents odor  
and propagation  
of bacteria



Compressive  
undersock provides  
edema control  
in the foot



<sup>1</sup>Inelastic compression legging produces gradient compression and significantly higher skin surface pressures compared with an elastic compression stocking. Kline CN, Macias BR, Kraus E, Neuschwander TB, Angle N, Bergan J, Hargens AR. Published in Vascular, Volume 16, No. 1, 2008

\*EZ-on panel is optional use.

# compression options for every patient



## Juxtalite indications

- Mild edema
- Active and healed venous ulcers
- Chronic Venous Insufficiency



## Juxtalite HD indications

- Moderate edema
- Active and healed venous ulcers
- Phlebolympheidema
- Chronic Venous insufficiency

**COMING SOON**  
black color option



**NEW EZ-on panel\***  
included with every  
Juxtalite HD




## Each Juxtalite kit comes with:

1 adjustable compression leg wrap, 2 compressive undersocks, 1 BPS card

circumference in centimeters								
size	S	M	M - FC	L	L - FC	XL	XL - FC	XXL
cC	26-36	33-44	44-54	41-51	51-61	37-48	48-61	53-64
cB	19-22	22-27	22-27	25-32	25-32	30-37	30-37	35-42
actual product length in centimeters								
size	short				long			
ℓB-D	28				33			
sizing measurement	short				long			
ℓ floor-D1	<44				>44			

Product lengths may differ from body measurements.



## mediven® dual layer compression system

- Elastic, two-layer stocking system
- System includes 2 closed-toe compressive liners and 2 open-toe compressive stockings

## mediven dual layer indications

- Mild edema
- Chronic Venous insufficiency
- Active and healed venous ulcers

Includes Easy-On foot slipper donning aid.

circumference	S	M	L	XL	XXL
calf	11" - 15" 28-38 cm	12¼" - 16½" 31-42 cm	13½" - 18" 34-46 cm	14½" - 19½" 37-50 cm	18" - 22" 46-56 cm
ankle	7" - 8½" 18-21.5 cm	8½" - 9¾" 21.5-25 cm	9¾" - 11" 25-28 cm	11" - 12½" 28-32 cm	12½" - 14½" 32-37 cm
leg length	STANDARD: 15"-19" (38-47.5 cm)				



## mediven dual layer stocking system

Item Description	Available Colors		Category	Size
	BEIGE	WHITE**		
<b>CLOSED TOE</b>				
D040 Liner Refill pair	–	•	1	S-XXL
D240 Calf 30-40mmhg* pairs	•	–	1	S-XXL
D240 Calf 30-40mmhg single - 10 pack	•	–	9	S-XXL
D340 Calf 40-50mmhg* pairs	•	–	1	S-XXL
D340 Calf 40-50mmhg singles	•	–	2	S-XXL

**Relevant Codes: HCPCS A6531** gradient compression stocking, below knee, 30-40 mmHg, each. Report with correct modifiers, RT (right) and/or LT (left). **HCPCS A6532** gradient compression stocking, below knee, 40-50 mmHg, each. Report with correct modifiers, RT (right) and/or LT (left). Code usage should be verified by responsible parties in accordance with all published guidelines.

\*Includes 2 liners and 2 outer stockings

\*\*white color code is 0 only for mediven dual layer liner refills

### circaid juxtalite lower leg

### circaid juxtalite hd with EZ-on panel lower leg

Beige	Compression wraps length (Short)
CJL1S001	small
CJL1S002	medium
CJL2S002	medium full calf
CJL1S003	large
CJL2S003	large full calf
CJL1S004	extra-large
CJL2S004	extra-large full calf
CJL1S005	XXL

Beige	Black	Compression wraps length (Short)
CJL3S001	CJL3S501	small
CJL3S002	CJL3S502	medium
CJL4S002	CJL4S502	medium full calf
CJL3S003	CJL3S503	large
CJL4S003	CJL4S503	large full calf
CJL3S004	CJL3S504	extra-large
CJL4S004	CJL4S504	extra-large full calf
CJL3S005	CJL3S505	XXL

Beige	Compression wraps length (Long)
CJL1L001	small
CJL1L002	medium
CJL2L002	medium full calf
CJL1L003	large
CJL2L003	large full calf
CJL1L004	extra-large
CJL2L004	extra-large full calf
CJL1L005	XXL

Beige	Black	Compression wraps length (Long)
CJL3L001	CJL3L501	small
CJL3L002	CJL3L502	medium
CJL4L002	CJL4L502	medium full calf
CJL3L003	CJL3L503	large
CJL4L003	CJL4L503	large full calf
CJL3L004	CJL3L504	extra-large
CJL4L004	CJL4L504	extra-large full calf
CJL3L005	CJL3L505	XXL



**System includes:**  
**1 juxtalite Legging,**  
**2 compressive undersocks**  
**and 1 BPS card.**

**System includes: 1 juxtalite hd Legging,**  
**2 compressive undersocks and 1 BPS card.**

**Relevant Codes: HCPCS A6545** for an open venous stasis ulcer, billed with the AW modifier. Limited to one compression wrap per six months, per leg with RT (right) and/or LT (left) modifiers. Gradient compression, non-elastic, below knee 30-50 mmHg. Code usage should be verified by responsible parties in accordance with all published guidelines.

### circaid compressive undersocks

Item	Description	Colors
		BLACK
CU001511	circaid compressive undersock 15-25 mmHg (standard)	•
CU001512	circaid compressive undersock 15-25 mmHg (large)	•
CU002511	circaid compressive undersock 25-35 mmHg (standard)	•
CU002512	circaid compressive undersock 25-35 mmHg (large)	•

Compressive undersock provides compression only to the foot, for greater edema control.

**medi USA**  
**6481 Franz Warner Pkwy**  
**Whitsett, NC 27377**  
**T 800-633-6334**  
**F 888-570-4554**  
**info@mediusa.com**  
**mediusa.com**